

1 **SUPPLEMENTARY METHODS**

2 **Tissue digestion and cell isolation**

3 Tumour tissues, tumour-associated lymph nodes, and colorectal healthy mucosa were
4 collected in IMDM+Glutamax medium (Gibco) complemented with 20% fetal calf
5 serum (FCS) (Greiner Bio-One) after surgical resection. Samples were cut into small
6 fragments in a petri dish, and enzymatically digested with 1 mg/mL collagenase D
7 (Roche Diagnostics) and 50 µg/mL DNase I (Roche Diagnostics) in 5 mL of
8 IMDM+Glutamax medium for 30 min at 37°C in gentleMACS C tubes (Miltenyi Biotec).
9 During and after incubation, cell suspensions were mechanically dissociated on the
10 gentleMACS Dissociator (Miltenyi Biotec). Cell suspensions were filtered through a
11 70-µm cell strainer (Corning) and washed in IMDM+Glutamax medium. Cell number
12 and viability was determined with the Muse Count & Viability Kit (Merck) on the Muse
13 Cell Analyser (Merck). Cells were cryopreserved based on the number of viable cells
14 in liquid nitrogen until time of analysis in 50% FCS and 10% dimethyl sulfoxide (DMSO)
15 (Merck). Peripheral blood mononuclear cells (PBMCs) were isolated by Ficoll-Paque
16 (provided by apothecary LUMC) density-gradient centrifugation and cryopreserved in
17 liquid nitrogen until time of analysis in 20% FCS and 10% DMSO.

18

19 **Mass cytometry antibodies**

20 Heavy metal isotope-tagged monoclonal antibodies are listed in [Table S2](#). Purified
21 antibodies were conjugated with heavy metal reporters in-house using the MaxPar X8
22 Antibody Labeling Kit (Fluidigm) according to the manufacturer's instructions. All
23 antibodies were titrated to determine the optimal labelling concentration.

24

25 **Mass cytometry antibody staining and data acquisition**

26 Percoll (GE Healthcare) density-gradient centrifugation was performed to isolate
27 immune cells from tumour tissues and colorectal healthy mucosa. Cells were washed
28 in Maxpar Cell Staining Buffer (CSB, Fluidigm), and incubated with 1 mL CSB
29 containing 1 μ M Cell-ID intercalator-¹⁰³Rh (Fluidigm) for 15 min at room temperature
30 (rT). Cells were washed, incubated with human Fc receptor block (BioLegend) for 10
31 min at rT, and stained with cell surface antibodies for 45 min at rT in a final volume of
32 100 μ L. After washing, cells were incubated with 1 mL Maxpar Fix and Perm buffer
33 (Fluidigm Sciences) containing 0.125 μ M Cell-ID intercalator-Ir (Fluidigm) overnight at
34 4°C. Cells were acquired on a Helios-upgraded CyTOF2 and Helios mass cytometer
35 (Fluidigm) at an event rate of <500 events/sec in de-ionized water containing 10x
36 diluted EQ Four Element Calibration Beads (Fluidigm). Data were normalized with the
37 normalization passport EQ-P13H2302_ver2 for each experiment.

38

39 **Mass cytometry data analysis**

40 Correlation analysis of the presence of immune cell clusters across samples (Figure
41 7) was performed using Spearman's rank correlation in the Corrplot R package. PD-1
42 subsets in CRCs (Figure S9) were gated based on a healthy control PBMC sample
43 into PD-1 negative cells (no PD-1 expression), PD-1 intermediate cells (PD-1
44 expression levels to a similar extent as healthy control PBMCs), and PD-1 high cells
45 (PD-1 expression levels surpassing that of healthy control PBMCs) as described in a
46 previous study on functional diversity between PD-1 subsets.¹ Tumours with less than
47 100 cells in one of the PD-1 subsets were excluded from this analysis. Sample and
48 immune cell cluster t-distributed Stochastic Neighbour Embedding (t-SNE) maps
49 (Figure S8) were computed as described previously.²

50

51 **Flow cytometry antibody staining**

52 Cells were first incubated with human Fc receptor block and cell surface antibodies for
53 45 min at 4°C. After washing, cells were stained for intracellular proteases and
54 cytokines using Fixation Buffer and Intracellular Staining Perm Wash Buffer
55 (BioLegend) or for FOXP3 expression using FOXP3 Transcription Factor Staining
56 Buffer Set (eBioscience). Compensation was carried out with CompBeads (BD
57 Biosciences) and ArC reactive beads (Life Technologies). Cells were acquired on a
58 LSR II and LSR Fortessa flow cytometer (BD Biosciences) running FACSDiva
59 software version 8.0 (BD Biosciences). Data were analysed with FlowJo software
60 version 10.2 (Tree Star Inc). Antibody details are listed in [Table S3](#).

61

62 **Immunohistochemical staining**

63 FFPE blocks from colorectal cancers (CRCs) were obtained from the department of
64 Pathology at the Leiden University Medical Centre (Leiden, The Netherlands). Tumour
65 mismatch repair (MMR) status (MMR-proficient or MMR-deficient) was determined by
66 immunohistochemical detection of PMS2 (anti-PMS2 antibodies; clone EP51, DAKO)
67 and MSH6 (anti-MSH6 antibodies; clone EPR3945, Abcam) proteins.³ MMR-
68 deficiency was determined by lack of expression of at least one of the MMR-proteins
69 in the presence of an internal positive control. Immunohistochemical detection of
70 human leukocyte antigen (HLA) class I expression was performed with HCA2 and
71 HC10 monoclonal antibodies (Nordic-MUbio), and classified as HLA class I positive,
72 weak, or loss as described previously⁴.

73

74 **Multispectral immunofluorescence detection**

75 Five- μ m frozen tissue sections were cut on adhesive immunohistochemistry slides.
76 Following fixation in ice-cold 100% methanol for 5 min and wash in PBS, the tissues
77 were incubated with Superblock buffer (Thermo Fisher Scientific) for 30 min at rT. For
78 the detection of CD127, TSA signal amplification was performed with Opal520 from
79 the Opal 7-colour manual IHC kit (Perkin Elmer), according to manufacturer's
80 instructions. Thereafter, the tissues were incubated with primary TCR V beta F1 and
81 CD45RO antibodies to be detected indirectly overnight at rT, followed by incubation
82 with corresponding fluorescent secondary antibodies CF680 and CF633, respectively,
83 for 1h at rT. The tissues were then incubated with directly conjugated primary CD3-
84 AF594 and CD7-AF647 antibodies for 5h at rT. CD3 was conjugated to AF594 using
85 the Alexa Fluor 594 antibody labelling kit (Thermo Fisher Scientific). Lastly, the tissues
86 were incubated with 1 μ M DAPI for nuclear counterstaining, and mounted with
87 Prolong® Gold Antifade Reagent (Cell Signaling Technology). Antibody details are
88 listed in [Table S4](#).

89

90 **References**

- 91 1. Thommen DS, Koelzer VH, Herzig P, et al. A transcriptionally and functionally
92 distinct PD-1(+) CD8(+) T cell pool with predictive potential in non-small-cell
93 lung cancer treated with PD-1 blockade. *Nat Med* 2018 doi: 10.1038/s41591-
94 018-0057-z
- 95 2. van Unen V, Li N, Molendijk I, et al. Mass Cytometry of the Human Mucosal
96 Immune System Identifies Tissue- and Disease-Associated Immune Subsets.
97 *Immunity* 2016;44(5):1227-39. doi: 10.1016/j.immuni.2016.04.014
- 98 3. Hall G, Clarkson A, Shi A, et al. Immunohistochemistry for PMS2 and MSH6 alone
99 can replace a four antibody panel for mismatch repair deficiency screening in
100 colorectal adenocarcinoma. *Pathology* 2010;42(5):409-13. doi:
101 10.3109/00313025.2010.493871 [published Online First: 2010/07/17]
- 102 4. Ijsselsteijn ME, Petitprez F, Lacroix L, et al. Revisiting immune escape in
103 colorectal cancer in the era of immunotherapy. *Br J Cancer* 2019;120(8):815-
104 18. doi: 10.1038/s41416-019-0421-x [published Online First: 2019/03/14]

105